
# Kobold®

## Kobold Prostate HDR Template Kobold Prostate HDR Stepper Template **INSTRUCTIONS FOR USE**


United States Patent No. D755,385  
www.koboldmedical.com

### Prostate HDR Templates

PART NUMBER	DESCRIPTION	QTY.
801-000	KOBOLD PROSTATE HDR TEMPLATE - 15 ga	1
802-000	KOBOLD PROSTATE HDR TEMPLATE - 16 ga/5 Fr	1
803-000	KOBOLD PROSTATE HDR TEMPLATE - 17 ga/4 Fr	1
804-000	KOBOLD PROSTATE HDR TEMPLATE - 18 ga	1
805-000	KOBOLD PROSTATE HDR TEMPLATE - 6 Fr	1


\*PICTURED ABOVE: 801-000 Kobold Prostate HDR Template - 15 ga

### Prostate HDR Stepper Templates

PART NUMBER	DESCRIPTION	QTY.
810-000	KOBOLD PROSTATE HDR STEPPER TEMPLATE - 15 ga	1
811-000	KOBOLD PROSTATE HDR STEPPER TEMPLATE - 16 ga/5 Fr	1
812-000	KOBOLD PROSTATE HDR STEPPER TEMPLATE - 17 ga/4 Fr	1
813-000	KOBOLD PROSTATE HDR STEPPER TEMPLATE - 6 Fr	1

\*PICTURED ABOVE: 810-000 Kobold Prostate HDR Stepper Template - 15 ga



Authorized Representative:  
Emergo Europe  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands

Manufactured For:  
Kobold, LLC  
23403 E. Mission Ave., Suite 220E  
Liberty Lake, WA 99019  
USA

# SYMBOLS GLOSSARY

	STD/REG REF	STANDARD/REGULATION TITLE	SYMBOL TITLE	SYMBOL EXPLANATION
	21 CFR 801.15(c) (1)(i)F  21 CFR 801.109	Labeling-Medical devices; prominence of required label statements.  Labeling-Prescription devices.	Prescription only	Caution: Federal law (U.S.) restricts this device to sale by or on the order of a physician
	ISO 15223-1, Clause 5.4.4	Medical devices. Symbols to be used with medical device labels, labelling and information to be supplied.	Attention: Read all warnings and precautions in instructions for use	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
	ISO 15223-1, Clause 5.2.7	Medical devices. Symbols to be used with medical device labels, labelling and information to be supplied.	Non sterile	Indicates a medical device that has not been subjected to a sterilization process.
	ISO 15223-1, Clause 5.2.8	Medical devices. Symbols to be used with medical device labels, labelling and information to be supplied.	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened.
	ISO 15223-1, Clause 5.4.2	Medical devices. Symbols to be used with medical device labels, labelling and information to be supplied.	Do not reuse	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	ASTM F2503	Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment.	Magnetic Resonance (MR) unsafe	Keep away from magnetic resonance imaging (MRI) equipment.
	ISO 15223-1, Clause 5.1.6	Medical devices. Symbols to be used with medical device labels, labelling and information to be supplied.	Catalogue or model number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	ISO 15223-1, Clause 5.1.5	Medical devices. Symbols to be used with medical device labels, labelling and information to be supplied.	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	ISO 15223-1, Clause 5.1.1	Medical devices. Symbols to be used with medical device labels, labelling and information to be supplied.	Manufacturer	Indicates the medical device manufacturer.
	765/2008/EC 768/2008/EC MDD 93/42/EEC Articles 4,11,12,17, Annex II	The requirements for accreditation and market surveillance relating to the marketing of products; Medical Device Directive.	CE marking	Signifies European technical conformity.
	ISO 15223-1, Clause 5.1.2	Medical devices. Symbols to be used with medical device labels, labelling and information to be supplied.	Authorized representative in the European Community	Indicates the authorized representative in the European Community.

**ENGLISH****DESCRIPTION**

The Kobold Prostate HDR Template and Stepper Template are US patented accessories to be used where ultrasound-guided interstitial temporary high dose rate (HDR) brachytherapy treatment of the prostate is accepted by up-to-date clinical guidelines. The templates are best used with biplanar transrectal ultrasound and confirmatory fluoroscopy.

The templates are single-use, disposable device shipped nonsterile. These devices are for continuous use of up to 24 hours of contact with the patient.

For workflows that require a stepper, the Kobold Prostate HDR Stepper Template fits into the Elekta Martinez Template Holder. The Prostate HDR template is ergonomically designed for freehand HDR technique.

It is customizable for almost any prostate size or shape. Large and small prostate sizes can all be accommodated with one template. The templates include four suture holes at the corners to secure the template to the perineum. It features rounded edges for patient comfort and reduction of pressure points on the skin.

This highly customizable template can be utilized for various clinical scenarios. It provides urethral sparing as well as the ability to implant suspected extracapsular disease and the proximal seminal vesicles, should that be desired. The ridged grip on the periphery of the template allows for easy handling during the procedure. Extra needle holes are provided to allow optimal implantation of irregular or asymmetric prostate glands. The templates are available in various gauge sizes to match the most common needles sizes. Easily identifiable grooves are marked on the template to indicate typical potential needle placement locations for most prostate sizes, and also help guide the needle into the hole.

**INDICATIONS FOR USE**

The Kobold Prostate HDR Template and Stepper Template are indicated for use as an accessory for high dose rate brachytherapy treatment of the prostate.

**CONTRAINDICATIONS**

As per clinical guidelines and published and accepted standard clinical practice.

**CAUTIONS**

Each physician is responsible for using the appropriate patient selection and the correct use of this template based on his/her own experience, training, and medical judgment. The physician must be trained in the proper use of the afterloader.

The Kobold Prostate HDR Template and the Kobold Prostate HDR Stepper Template have not been tested for MR compatibility and safety.

**WARNINGS**

Possible complications associated with the use of this device are related to needle insertion and usage techniques of the responsible physician.

The template must be properly sterilized prior to use.

Using a damaged template can lead to incorrect treatment of the patient and may cause radiation hazards to patients and medical staff. Physician users should familiarize themselves with the feel of the needle and template fit prior to any procedure as there may be differences in materials or slight gauge differences between various needle manufacturers.

**DELIVERY CONDITION AND TESTING**

Each Kobold Prostate HDR Template and HDR Stepper Template is carefully packaged. It is supplied in a NONSTERILE condition. The packaging should be unopened and undamaged upon receipt. Remove the components from the packaging in a clean environment. Prior to use, the user must ensure that the template is undamaged.

**PARTS/CONDITION OF PARTS**

If the proper condition of this product or part is in question, it must not be utilized, and the part should be returned immediately for inspection by the manufacturer. If damage to the template is obvious, please contact the manufacturer for replacement part.

The manufacturer is not responsible for the use of damaged or modified parts/products or the use of products that do not meet factory specifications.

**STERILIZATION**

The hospital or clinic is responsible for in-house procedures for the inspection and packaging of the devices in a manner that will ensure adequate steam sterilant penetration and adequate drying.

Double pouch in sterilization pouch and seal.

See the table below for recommended steam sterilization parameters that have been validated by Kobold to provide a 10<sup>-6</sup> sterility assurance level (SAL).

Table 1 Standard Autoclave Specifications

Cycle Type	Temperature	Exposure Time	Dry Time
Pre-vacuum	132°C 270°F	4 minutes	20 minutes

Note: The sterilizer manufacturer's instructions for operation and load configuration should be followed explicitly.

**SETUP AND USE**

The Kobold Prostate HDR Template and HDR Stepper Template are best used with a biplanar transrectal ultrasound probe and confirmatory fluoroscopy. The type of anesthesia and operating room setting is determined by the individual physician. Standard sterile technique is utilized throughout the procedure.

- Before using this single-use, disposable template, follow the Sterilization Instructions.
- If choosing to perform the procedure with a Stepper, the Kobold Prostate HDR Stepper Template is designed to be utilized with the Elekta Martinez Template Holder. The template is placed into the holder with the raised "Kobold" label (front) facing the physician. The side of the template with deep grooves (back, image B), faces the patient. The two locking pins are slid down the channels securing the template into place as seen in images A and B below.

Note: The grooves and wells on the back of the template are only for manufacturing purposes and have no clinical function.

- A radioopaque Foley catheter is placed in the urethra/bladder (under standard sterile conditions). Standard transrectal ultrasound measurements of the prostate are accomplished. It is recommended that the penis and scrotum are held up out of the field with an adhesive dressing. Typically, four gold seed fiducial markers are placed in the prostate, two at the apex, two at the base, or per physician preference.
- The physician presses on the perineum and feels the edge of the symphysis pubis. This is used to estimate placement of the template on the perineum, ensuring there is adequate space for the ultrasound probe posteriorly. The template is notched to fit any size ultrasound probe.
- The templates are available in various Gauge/French sizes (see #7). The 16 gauge template typically is to be used with 16 gauge needles, and the 15 gauge template is to be used with 15 gauge needles and so on. It is up to individual physician preference to use plastic or metal needles, and these should be used according to their respective instructions for use.

Caution: Depending on the needle manufacture, the needles may be smaller or larger than standard gauge and French sizes.

- The holes for the needles are intended to be relatively tight with a dry needle to limit slippage after the implant is complete. IMPORTANT: DIPPING THE TIP OF THE NEEDLE IN ISOPROPYL ALCOHOL PRIOR TO INSERTION will allow the needle to slide easily through the hole. When the alcohol evaporates, the needle should be relatively secure. It is important not to contaminate the needle holes with lubricating jelly or other fluids as this may cause the needles to slide too easily in the holes. Make certain that no fluid enters the lumen of the needle.

- Some physicians prefer that the needles are held very tightly by the template after the implant is complete. If this is desired, then using a template with a gauge size one step smaller than the needle gauge size may provide the desired post-implant stability of the needles within the template (for example 16 G needle may be held tighter in a 17 G template). In this case, IT IS IMPERATIVE TO DIP THE TIP OF THE NEEDLE IN ISOPROPYL ALCOHOL PRIOR TO INSERTION.

Note: It is imperative to test the needle/template combination in a simulated environment prior to performing an implant on a patient to assure desired fit of the needles in the template and the relative ease or difficulty of adjustment of depth if required. Again, there may be differences in materials or slight gauge differences between various needle manufacturers.

- The template has a circular groove that connects a ring of needle holes. This guide corresponds to the typical positioning of the peripheral needles, implanted at the prostate capsule, for most cases (i.e., prostate sizes of 25 to 50 cc). Larger glands or asymmetric glands may need to use the outer ring of needle holes. For smaller glands one would use the inner concentric rings.

The inner "U" shaped groove is also a guide for potential needle locations of the inner needles for most prostate sizes. Usually two essentially concentric rings of needles are adequate for most glands. This allows excellent sparing of the urethra, yet optimal peripheral dose and improved dose homogeneity.

Typical implants require 16 to 19 needles. This is highly dependent on physician preference and technique, and the template is customizable for most shapes and sizes of prostate glands and individual physician preference. See figures C, D, E, and F for needle locations for examples of typical implants.

Both types of templates have the exact same needle hole positions.

9. The physician places peripheral needles carefully in the prostate gland using transrectal ultrasound to verify needle depth covers the entire base of the gland. The anterior and peripheral needles are placed first so as to not obscure ultrasound visualization for the subsequent needle placement.

10. The central needles are placed next, usually starting at the anterior prostate and working posteriorly.

11. Finally the most posterior needles are placed. According to the clinical scenario and physician preference, the posterior needles may be used to implant the proximal seminal vesicles.

12. Care is taken to avoid the rectum and bladder. In the sagittal plane, the physician should verify the depth of the needles in conjunction with fluoroscopy and fiducial markers to ensure coverage of the base of the prostate gland according to physician discretion and clinical scenario.

13. When all needles are placed, the template is equipped with four suture holes for securing the template to the perineum (with or without bolster/padding according to physician preference).

14. The template/needles are numbered according to institution preference. Each needle is assigned a number and these are recorded for future reference. A sample template numbering guide can be found at [www.koboldmedical.com](http://www.koboldmedical.com) and may be printed for each patient's implant. It is imperative that the numbering system of the catheters is consistent throughout all phases of the procedure. This should be re-verified by at least 2 independent, qualified clinical personnel prior to treatment.

15. The physician may choose to record the measurement from the surface of the template to the end of each needle for another reference verification measurement together with fiducial/needle relationship to help ensure proper needle positioning. If post-implant adjustment of a catheter is needed, slightly twisting the needle should help release the template grip.

Tip: One may use a permanent ink surgical marker and mark the needle at the entrance point into the template for easy reference of potential needle shifting after implant and prior to treatment.

16. Real-time ultrasound treatment planning or CT simulation and treatment planning according to institution and physician preference then follows. The number of fractions, catheter position and numbering verification, and other quality assurance procedures are determined by the treating physician, institution, and standard practices.

Note: It is imperative that each physician verifies the needle position and depth (coverage of the base of the prostate) during simulation, planning, and treatment. One must ensure that the needles have not shifted position between phases of the treatment. Needle adjustments prior to treatment may need to be made.

Note: If parts supplied by other manufacturers are utilized in combination with Kobold products, the physician is responsible for incompatibilities and is responsible for following other manufacturers' instructions for use.

## STORAGE AND DISPOSAL

Before use, store the template in the Kobold-provided wrapping/bag. Store at an ambient temperature (15° to 30°C). Safely dispose template after single use according to standard practice.

## WARRANTY

The Kobold Prostate HDR Template and Kobold Prostate HDR Stepper Template are warranted to be free from defects in workmanship and materials, under normal use, for a period of 30 days from the date of original purchase from Kobold. The sole obligation of Kobold under this warranty shall be, at the company's option, to replace the product or part with a comparable product or part.

Warranty is not valid if part is damaged due to customer's mishandling or misuse of said part.

Kobold makes no warranty, whatsoever, to any persons other than the original customer.

## KOBOLD TECHNICAL SUPPORT

For technical questions relating to this product, its use, or for reporting any incidents involving the product, contact Kobold Technical Support:  
+1 509-703-5090 | www.koboldmedical.com | info@koboldmedical.com

Image A

Image B

Image C - Axial CT image of typical implant needle locations with 17 metal needles.


Image D - Axial CT image of typical implant needle locations with 19 plastic needles.

Image E - Prostate HDR Template with typical implant needle locations with 17 needles.

Image F - Prostate HDR Stepper Template with typical implant needle locations with 19 needles.

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