# Kobold®

Kobold Henschke Tandem and Ovoid, Tandem and Ring, Fletcher Tandem and Ovoid,  
and Vaginal Cylinder Applicator Sets

**REF** 704, 707, 708, 715 Series

## INSTRUCTIONS FOR USE

[www.koboldmedical.com](http://www.koboldmedical.com)

Manufactured For:  
Kobold, LLC  
23403 E. Mission Ave., Suite 220E  
Liberty Lake, WA 99019  
USA

Authorized Representative:  
Emergo Europe  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands

# TABLE OF CONTENTS


|                                                           |          |
|-----------------------------------------------------------|----------|
| <b>Symbols Glossary.....</b>                              | <b>i</b> |
| <b>Applicator Sets – Contents.....</b>                    | <b>1</b> |
| Kobold® Henschke-model Tandem and Ovoid Set.....          | 1        |
| Kobold® Tandem and Ring Applicator Set (series).....      | 2        |
| Kobold® Fletcher-model Tandem and Ovoid Set (series)..... | 3        |
| Kobold® Vaginal Cylinder Applicator Set.....              | 3        |
| <b>Instructions for Use.....</b>                          | <b>4</b> |
| English.....                                              | 4        |
| Deutsch.....                                              | 6        |
| Français.....                                             | 9        |
| Español.....                                              | 12       |
| Italiano.....                                             | 15       |
| Português.....                                            | 18       |
| Ελληνικά.....                                             | 20       |
| Svenska.....                                              | 24       |
| Norsk.....                                                | 26       |
| Dansk.....                                                | 29       |
| Suomi.....                                                | 31       |
| Nederlands.....                                           | 34       |

# SYMBOLS GLOSSARY

| STD/REG REF                                                                                                                                                      | STANDARD/REGULATION TITLE                                                                                                   | SYMBOL TITLE                                                         | SYMBOL EXPLANATION                                                                                                                                                                                                          |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <br>21 CFR 801.15(c)(1)(i)F<br>21 CFR 801.109                                   | Labeling-Medical devices; prominence of required label statements.<br><br>Labeling-Prescription devices.                    | Prescription only                                                    | Caution: Federal law (U.S.) restricts this device to sale by or on the order of a physician                                                                                                                                 |
| <br>ISO 15223-1, Clause 5.4.4                                                   | Medical devices. Symbols to be used with medical device labels, labelling and information to be supplied.                   | Attention: Read all warnings and precautions in instructions for use | Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself. |
| <br>ISO 15223-1, Clause 5.2.7                                                   | Medical devices. Symbols to be used with medical device labels, labelling and information to be supplied.                   | Non-sterile                                                          | Indicates a medical device that has not been subjected to a sterilization process.                                                                                                                                          |
| <br>ISO 15223-1, Clause 5.2.8                                                   | Medical devices. Symbols to be used with medical device labels, labelling and information to be supplied.                   | Do not use if package is damaged                                     | Indicates a medical device that should not be used if the package has been damaged or opened.                                                                                                                               |
| <br>ISO 15223-1, Clause 5.4.2                                                   | Medical devices. Symbols to be used with medical device labels, labelling and information to be supplied.                   | Do not reuse                                                         | Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.                                                                                                          |
| <br>ASTM F2503                                                                 | Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment.             | Magnetic Resonance (MR) unsafe                                       | Keep away from magnetic resonance imaging (MRI) equipment.                                                                                                                                                                  |
| <br>ISO 15223-1, Clause 5.1.6                                                 | Medical devices. Symbols to be used with medical device labels, labelling and information to be supplied.                   | Catalogue or model number                                            | Indicates the manufacturer's catalogue number so that the medical device can be identified.                                                                                                                                 |
| <br>ISO 15223-1, Clause 5.1.5                                                 | Medical devices. Symbols to be used with medical device labels, labelling and information to be supplied.                   | Batch code                                                           | Indicates the manufacturer's batch code so that the batch or lot can be identified.                                                                                                                                         |
| <br>ISO 15223-1, Clause 5.1.1                                                 | Medical devices. Symbols to be used with medical device labels, labelling and information to be supplied.                   | Manufacturer                                                         | Indicates the medical device manufacturer.                                                                                                                                                                                  |
| <br>765/2008/EC<br>768/2008/EC MDD<br>93/42/EEC Articles 4,11,12,17, Annex II | The requirements for accreditation and market surveillance relating to the marketing of products; Medical Device Directive. | CE Marking<br>Conformité Européenne Notified Body Reference no. #### | Signifies European technical conformity.                                                                                                                                                                                    |
| <br>ISO 15223-1, Clause 5.1.2                                                 | Medical devices. Symbols to be used with medical device labels, labelling and information to be supplied.                   | Authorized representative in the European Community                  | Indicates the authorized representative in the European Community.                                                                                                                                                          |

# KOBOLD APPLICATOR SETS - CONTENTS

Note: Some pictured parts are not to scale. If included in Set, Mini Stainless Steel Hex Driver and Stainless Steel Ruler are not pictured.


## 704-000 Henschke Tandem and Ovoid Applicator Set

| ITEM NO. | PART NUMBER | DESCRIPTION                           | QTY. |
|----------|-------------|---------------------------------------|------|
| 1        | 704-101     | HENSCHKE OVOID, 2.0CM DIAMETER        | 2    |
| 2        | 704-102     | HENSCHKE OVOID, 2.5CM DIAMETER        | 2    |
| 3        | 704-103     | HENSCHKE OVOID, 3.0CM DIAMETER        | 2    |
| 4        | 704-106     | HENSCHKE COLPOSTAT - PATIENT LEFT     | 1    |
| 5        | 704-104     | HENSCHKE COLPOSTAT - PATIENT RIGHT    | 1    |
| 6        | 602-105     | TANDEM STRAIGHT                       | 1    |
| 7        | 704-109     | HENSCHKE TANDEM - CURVE 25 DEGREE     | 1    |
| 8        | 704-111     | HENSCHKE TANDEM - CURVE 40 DEGREE     | 1    |
| 9        | 704-113     | HENSCHKE TANDEM - CURVE 60 DEGREE     | 1    |
| 10       | 704-114     | HENSCHKE TANDEM CLAMP                 | 1    |
| 11       | 601-252     | CERVICAL STOP                         | 1    |
| 12       | 602-108     | SPACING BRACKET PATIENT RIGHT         | 1    |
| 13       | 602-109     | SPACING BRACKET PATIENT LEFT          | 1    |
| 14       | 601-256     | 4MM LOCK KNOB                         | 2    |
| 15       | 601-255     | HEX SCREW M3 X 0.5MM                  | 14   |
| 16       | 602-111     | STERILIZATION CAP O.D. 0.25           | 30   |
| 17       | 704-118     | HENSCHKE-MODEL SET STERILIZATION TRAY | 1    |
| 18       | 601-266     | SCREW BOX                             | 1    |
| 19       | 601-262     | MINI STAINLESS STEEL HEX DRIVER       | 1    |
| 20       | 601-263     | STAINLESS STEEL RULER                 | 1    |

United States Patent No. 9,180,311  
United States Patent No. 9,895,551


707-001 Tandem and 30°/45°/60° Ring Applicator Set

| ITEM NO. | PART NUMBER | DESCRIPTION                            | QTY. |
|----------|-------------|----------------------------------------|------|
| 1        | 707-101     | 30 DEGREE X 6 CM TANDEM                | 1    |
| 2        | 707-103     | 45 DEGREE X 6 CM TANDEM                | 1    |
| 3        | 707-105     | 60 DEGREE X 6 CM TANDEM                | 1    |
| 4        | 707-107     | 30 DEGREE X 4 CM TANDEM                | 1    |
| 5        | 707-109     | 45 DEGREE X 4 CM TANDEM                | 1    |
| 6        | 707-111     | 60 DEGREE X 4 CM TANDEM                | 1    |
| 7        | 707-113     | 30 DEGREE X 2 CM TANDEM                | 1    |
| 8        | 707-115     | 45 DEGREE X 2 CM TANDEM                | 1    |
| 9        | 707-117     | 60 DEGREE X 2 CM TANDEM                | 1    |
| 10       | 707-119     | RING APPLICATOR, 30 DEG 28MM DIA       | 1    |
| 11       | 707-121     | RING APPLICATOR, 45 DEG 28MM DIA       | 1    |
| 12       | 707-122     | RING APPLICATOR, 60 DEG 28MM DIA       | 1    |
| 13       | 707-125     | RING BRACKET                           | 1    |
| 14       | 601-256     | 4MM LOCK KNOB                          | 2    |
| 15       | 707-126     | RING CAP 28MM - TOP                    | 1    |
| 16       | 707-127     | RING CAP 28MM - BOTTOM                 | 1    |
| 17       | 707-129     | RING APPLICATOR SET STERILIZATION TRAY | 1    |
| 18       | 602-111     | STERILIZATION CAP O.D. 0.25            | 30   |

Alternative Tandem and Ring Applicator Set Options

| PART NUMBER                   | DESCRIPTION                      |
|-------------------------------|----------------------------------|
| 707-002 TANDEM & 60° RING SET | ITEM NUMBERS: 1, 4, 7, 12, 13-18 |
| 707-003 TANDEM & 45° RING SET | ITEM NUMBERS: 2, 5, 8, 11, 13-18 |
| 707-004 TANDEM & 30° RING SET | ITEM NUMBERS: 3, 6, 9, 10, 13-18 |


## 708-001 Fletcher Tandem and Barrel Ovoid Applicator Set

| ITEM NO. | PART NUMBER | DESCRIPTION                           | QTY. |
|----------|-------------|---------------------------------------|------|
| 1        | 708-117     | FLETCHER BARREL OVOID, 1.5CM DIAMETER | 2    |
| 2        | 708-118     | FLETCHER BARREL OVOID, 2.0CM DIAMETER | 2    |
| 3        | 708-119     | FLETCHER BARREL OVOID, 2.5CM DIAMETER | 2    |
| 4        | 708-120     | FLETCHER BARREL OVOID, 3.0CM DIAMETER | 2    |
| 5        | 708-107     | FLETCHER COLPOSTAT - PATIENT LEFT     | 1    |
| 6        | 708-109     | FLETCHER COLPOSTAT - PATIENT RIGHT    | 1    |
| 7        | 602-105     | TANDEM STRAIGHT                       | 1    |
| 8        | 708-110     | FLETCHER TANDEM - CURVE 35 DEGREE     | 1    |
| 9        | 708-112     | FLETCHER TANDEM - CURVE 50 DEGREE     | 1    |
| 10       | 708-114     | FLETCHER TANDEM - CURVE 70 DEGREE     | 1    |
| 11       | 708-121     | FLETCHER TANDEM CLAMP                 | 1    |
| 12       | 601-252     | CERVICAL STOP                         | 1    |
| 13       | 602-108     | SPACING BRACKET PATIENT RIGHT         | 1    |
| 14       | 602-109     | SPACING BRACKET PATIENT LEFT          | 1    |
| 15       | 601-256     | 4MM LOCK KNOB                         | 2    |
| 16       | 601-255     | HEX SCREW M3 X 0.5MM                  | 8    |
| 17       | 602-111     | STERILIZATION CAP O.D. 0.25           | 30   |
| 18       | 708-116     | FLETCHER-MODEL SET STERILIZATION TRAY | 1    |
| 19       | 601-266     | SCREW BOX                             | 1    |
| 20       | 601-262     | MINI STAINLESS STEEL HEX DRIVER       | 1    |
| 21       | 601-263     | STAINLESS STEEL RULER                 | 1    |

United States Patent No. 9,180,311  
United States Patent No. 9,895,551


\*708-001 includes Item Numbers 1-21

## 708-002 Fletcher Tandem and Cylinder Ovoid Applicator Set

| ITEM NO. | PART NUMBER | DESCRIPTION                    | QTY. |
|----------|-------------|--------------------------------|------|
| 22       | 708-101     | FLETCHER OVOID, 1.0CM DIAMETER | 2    |
| 23       | 708-102     | FLETCHER OVOID, 1.5CM DIAMETER | 2    |
| 24       | 708-103     | FLETCHER OVOID, 2.0CM DIAMETER | 2    |
| 25       | 708-104     | FLETCHER OVOID, 2.5CM DIAMETER | 2    |
| 26       | 708-105     | FLETCHER OVOID, 3.0CM DIAMETER | 2    |

\*708-002 includes Item Numbers 5-26

## 715-000 Vaginal Cylinder Applicator Set

| ITEM NO. | PART NUMBER | DESCRIPTION                             | QTY. |
|----------|-------------|-----------------------------------------|------|
| 1        | 715-101     | VAGINAL CYLINDER OBTURATOR              | 1    |
| 2        | 715-103     | VAGINAL CYLINDER OBTURATOR COLLET       | 4    |
| 3        | 715-104     | VAGINAL CYLINDER OBTURATOR COLLET CAP   | 4    |
| 4        | 715-106     | 2.0CM DIAM. VAGINAL CYLINDER            | 1    |
| 5        | 715-108     | 2.5CM DIAM. VAGINAL CYLINDER            | 1    |
| 6        | 715-110     | 3.0CM DIAM. VAGINAL CYLINDER            | 1    |
| 7        | 715-112     | 3.5CM DIAM. VAGINAL CYLINDER            | 1    |
| 8        | 715-114     | VAGINAL CYLINDER SET STERILIZATION TRAY | 1    |
| 9        | 602-111     | STERILIZATION CAP O.D. 0.25             | 30   |

## INSTRUCTIONS FOR USE

## Description

*Kobold Henschke Tandem and Ovoid Applicator Set*  
*Kobold Tandem and Ring Applicator Set*  
*Kobold Fletcher Tandem and Ovoid Applicator Set*  
*Kobold Vaginal Cylinder Applicator Set*

The HDR **Kobold® Henschke-model Tandem and Ovoid Applicator Set** design is a modified Henschke type GYN applicator with a click fit connector interface. Kobold manufactures a version compatible with Elekta HDR afterloaders. This applicator is mainly utilized for HDR brachytherapy boosts to the cervix, uterus, and parametrial tissues in conjunction with external beam radiation therapy. This GYN set consists of precision manufactured intrauterine tandems and a pair of interlocking vaginal colpostats with Henschke type unshielded ovoids (all components are designed to allow for customization of the treated geometry). The distance between the ovoids is adjustable with a lock knob. Four differing tandem geometries are provided in the set: one straight, one with slight curvature (25°), one with medium curvature (40°), and one with maximal curvature (60°). These allow for maximal accommodation of the uterine shape and flexure.

The HDR **Kobold® Fletcher-model Tandem and Ovoid Applicator Set** is a modified Fletcher type GYN applicator with a click fit connector interface. Kobold manufactures versions compatible with Elekta HDR afterloaders. This applicator is mainly utilized for HDR brachytherapy boosts to the cervix, uterus, and parametrial tissues in conjunction with external beam radiation therapy. This GYN set consists of precision manufactured intrauterine tandems and a pair of interlocking vaginal colpostats with Fletcher type unshielded barrel or cylinder ovoids (all components are designed to allow for customization of the treated geometry). The distance between the ovoids is adjustable with a lock knob. Four differing tandem geometries are provided in the set: one straight, one with slight curvature (35°), one with medium curvature (50°), and one with maximal curvature (70°). These allow for maximal accommodation of the uterine shape and flexure.

The HDR **Kobold® Tandem and Ring Applicator Set** design is a modified Ring/Tandem type applicator with a click fit connector interface. Kobold manufactures versions compatible with Elekta HDR afterloaders. This applicator is particularly useful in providing an HDR brachytherapy boost to the cervix, uterus, and parametrial tissues in conjunction with external beam radiation therapy. The applicator consists of precision manufactured ring, ring cap, and intrauterine tandems, interlocking via a ring bracket. All components are designed for a defined geometry. The rings are provided in three different angles (30°, 45°, and 60°) with three different tandem end lengths (2 cm, 4 cm, and 6 cm) provided for each ring angle. One ring cap is provided.

The HDR **Kobold® Vaginal Cylinder Applicator Set** design is based on conventional cylindrical HDR delivery systems with a connector interface compatible with Elekta HDR afterloaders. This applicator is particularly useful in providing HDR radiation as definitive treatment, adjuvant treatment, or as a treatment boost to the vaginal cuff or vaginal walls, with or without external beam radiation therapy. The set consists of precision manufactured vaginal cylinders of varied diameters (2.0, 2.5, 3.0, 3.5 cm) for optimal radiation treatment delivery and anatomical customization.

The Applicator Sets are CT Compatible but MR Unsafe.

## Indications for Use

*Kobold Fletcher Tandem and Ovoid Applicator Set*  
*Kobold Henschke Tandem and Ovoid Applicator Set*  
*Kobold Tandem and Ring Applicator Set*  
*Kobold Vaginal Cylinder Applicator Set*

The CT Compatible/MR Unsafe HDR Kobold® Fletcher Tandem/Ovoid, Henschke Tandem/Ovoid, Tandem/Ring, and Vaginal Cylinder Applicator Sets are indicated for use in any case where high dose rate (HDR) brachytherapy treatment of the cervix/uterus (with the Tandem/Ovoid or Tandem/Ring type applicator) or the vagina (with the Vaginal Cylinder Applicator) are accepted by up-to-date clinical guidelines.

The Kobold® Applicator Sets are intended for continuous use of up to 24 hours of contact with patient.

## Contraindications

As per clinical guidelines and published and accepted standard clinical practice.

## Cautions

Each physician is responsible for using the appropriate technique and deciding on the applicator set to be used based on his/her own experience, training, and medical judgment. Each treatment application and prescribed radiation dose should be based on accepted published standard guidelines. The physician must be trained in the proper use of the Elekta afterloader.

Do not overtighten the hex screws (if applicable) as this may result in damage to the hex driver and may strip the threads of the titanium parts.

Ensure the sterilization caps are fitted to the tandem and colpostats, tandem and ring, or obturator at all times, except when the transfer tube is attached. This prevents solutions or foreign materials from entering the tandem/colpostats, tandem/ring, or obturator. Liquids, blood, or other kinds of contaminants entering the tubes can lead to malfunction of or damage to the Elekta afterloader and may obstruct the passage of the source wire.

The sterilization cap is a single use device that must be disposed of after each use. A reused cap might not securely close the tandem, colpostat, ring, or obturator so moisture could enter which could lead to malfunction or damage the Elekta afterloader.

Limitations of use are based on the maximal radiation dose acceptable to various organs and according to standard treatment protocols.

Due to the small diameter of the 2 cm vaginal cylinder and source distance from the tip of the source wire, there may be limitation on uniform prescription at the vaginal apex with the smallest cylinder (2 cm). Utilization of the largest vaginal cylinder that fits comfortably is recommended.

Incorrect length measurement can lead to patient misadministration.

Using a damaged part can lead to incorrect treatment of the patient and may cause radiation hazards to patients and medical staff.

When treating in a sterile environment (i.e., operating room), standard sterile technique is utilized

throughout the procedure.

Take care not to bend the tandems/obturator and colpostats/rings during use. The 4mm lock knobs are designed to be hand tightened. Do not utilize tools such as wrenches for tightening as this can damage the threading.

## Warnings

The Applicator Sets should only be used with the Elekta microSelectron® and Elekta Flexitron® afterloaders.

Inaccuracies in dwell positions within the tandem/obturator or in the colpostats/ovoids can lead to unintended dose delivery if the applicator is not characterized with proper physics commissioning and dose verification. Users must characterize the Applicator Sets through autoradiography to ensure accurate dwell positions when used with the Elekta afterloader.

The Applicator Sets must be properly cleaned and sterilized prior to each use.

Transfer tubes should be maintained free of twists, cracks, folds, or kinks and protected against the introduction of any fluid or foreign materials into the lumen.

Interchanging transfer tube connections will lead to incorrect treatment of the patient. Elekta applicators are equipped with matching connectors to prevent this from occurring.

Prior to each use, the physical condition of the titanium components should be free of cracks, kinks, folds, or corrosion. Plastic components should be free of cracks, discolorations, chalking, or folds. If any of these conditions exist, the device has reached the end of its useful life and must be replaced. Contact Kobold for replacement parts or servicing of parts.

The sterilization cap must remain in place until just prior to connection of the transfer tubes to the applicator, and a new one must be put in place immediately after treatment.

Possible complications associated with the use of this device are related to insertion and usage techniques and may include but are not limited to vaginal or uterine injury, vaginal foreshortening, chronic fibrosis, need for corrective surgery, chronic vaginal dryness and/or irritation, fistula, bowel or bladder injury, organ perforation, bleeding, infection, and second malignancy risk.

An applicator that is not completely disassembled might become damaged beyond repair during the cleaning and sterilization cycle or may not become fully sterile.

**This Instructions for Use document applies to the 704, 707, 708, 715 Applicator Set series only.**

**These Kobold Applicator Sets are compatible with the following Elekta afterloaders:**

**microSelectron® and Flexitron®.** Refer to the instructions for use for the appropriate Elekta transfer tube set when using this applicator set with the Elekta afterloader. Use the appropriate Elekta X-ray markers for taking radiographs, if applicable.

The applicator set should only be used by a physician who has received full training in the appropriate Elekta afterloader system. Use only Elekta transfer tubes to connect the applicator to the afterloader.

## Delivery Condition and Testing

Each Kobold Applicator Set is carefully packaged. It is supplied in a **nonsterile** condition. The packaging should be unopened and undamaged upon receipt. Remove the components from the packaging in a clean environment. These reusable applicators must be handled carefully to avoid damage. Prior to each treatment, the user must ensure that the applicator set is undamaged, contains all components (see **Applicator Sets – Contents** section) and still meets specifications. Each set must be properly commissioned by a physics team to verify applicator geometry and radiation dose delivery prior to clinical use.

Before each treatment, after connecting the transfer tube, use the applicable source position ruler for the afterloader to check the correct transfer tube/applicator combination length. For the Elekta microSelectron®, this test should be performed using the source position simulator.

## Cleaning and Sterilizing

**Note:** The entire Applicator Set must be fully disassembled before cleaning.

All parts of the Kobold Applicator Sets are reusable (except for the sterilization caps). For sterilization and use instructions for all Elekta products, refer to their corresponding instructions for use.

- These instructions apply only to products originally manufactured by Kobold, LLC.
- The user must clean and sterilize the applicator prior to the first and any additional use. It is critical that instruments are thoroughly manually cleaned before sterilization processing is performed.
- User defined cleaning and sterilization methods other than those described in these instructions must be validated by the user.
- Any applicator not listed above should be handled per the individual Instructions for Use.
- Ensure that the sterilization equipment meets the manufacturer's specifications.
- While cleaning, make sure that no liquid or other kind of contamination enters the applicators.
- The sterilization caps should be in place on the tandems/obturators and ring applicators/colpostats during cleaning and sterilization to prevent foreign materials from entering the tubes. Sterilization caps are single use parts.

## Cleaning

| STEP | CLEANING CRITERIA                                                                                                                                                                                                                      |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1    | Completely submerge instruments in neutral pH enzymatic cleaning solution, such as Enzol, at the concentration per manufacturer's instructions for use. Soak for a minimum of 20 minutes.                                              |
| 2    | If necessary use a soft-bristled, nylon brush to gently scrub the device until all visible soil has been removed. Particular attention must be given to the crevices, mated surfaces, connectors, and other hard to clean areas.       |
| 3    | Remove the device from the enzyme solution and rinse in flowing distilled water one (1) full minute for a minimum of three (3) times. Thoroughly and aggressively flush holes and other difficult-to-reach areas. Do not flush lumens. |
| 4    | If necessary, place prepared cleaning agents in a sonication unit. Completely submerge the device in a cleaning solution and sonicate for ten (10) minutes at 45-50 kHz.                                                               |
| 5    | Remove the device from the sonicator and rinse in flowing distilled water for one (1) full minute for a minimum of three (3) times. Thoroughly and aggressively flush holes and other difficult to reach areas (do not flush lumens).  |
| 6    | Inspect each unit for cleanliness and physical integrity. Repeat cleaning steps if additional cleaning is required.                                                                                                                    |

Steam Sterilization

- All parts of the applicator set can be steam sterilized.
- Sterilization caps should be in place on the tandems/obturator and ring applicators/colpostats during cleaning and sterilization to prevent foreign materials from entering the tubes.
- Remove the sterilization caps used during cleaning and place **new** sterilization caps on the tandems/obturator and ring applicators/colpostats **prior to sterilization**.
- Moist heat/steam is the preferred and recommended method for Kobold devices.
- See Table 1 for recommended steam sterilization parameters that have been validated by Kobold to provide a 10-6 sterility assurance level (SAL).
- The hospital/clinic is responsible for in-house procedures for reassembly, inspection, and packaging of the instruments after they are thoroughly cleaned in a manner that will ensure adequate steam sterilant penetration and adequate drying.
- Acceptable methods for wrapping instrument cases are described in EN 868 and AAMI ST79.

Table 1 Standard Autoclave Specifications

| CYCLE TYPE | TEMPERATURE | PRESSURE        | EXPOSURE TIME | DRY TIME   |
|------------|-------------|-----------------|---------------|------------|
| Pre-vacuum | 132°C       | 3 bar<br>44 psi | 4 minutes     | 15 minutes |
|            | 270°F       |                 |               |            |
| Pre-vacuum | 134°C       |                 | 3 minutes     | 20 minutes |
|            | 273°F       |                 |               |            |

- Sterilization manufacturer's recommendations should always be followed. When sterilizing multiple instruments sets in one sterilization cycle, ensure the manufacturer's maximum load is not exceeded.
- Applicator sets should be properly prepared and packaged in trays and/or cases that will allow steam to penetrate and make direct contact with all surfaces. Standard sterilization trays shall be double wrapped using no more than two layers of steam/hydrogen peroxide permeable hospital wrap. Rigid containers with vent systems require replacement of vent pads per the containers manufacturer's instructions.

Assembly and Use

*Kobold Fletcher Tandem and Ovoid Applicator Set*  
*Kobold Henschke Tandem and Ovoid Applicator Set*  
*Kobold Tandem and Ring Applicator Set*

Assembly Instructions

The sterilization caps should be placed on the attachment ends after assembly. It is recommended that the sterilization caps not be removed until the applicator is positioned in the patient and removal is necessary to connect to the HDR afterloader device.

For demonstration purposes only, the Kobold website provides links to videos of clinician assembling these applicator sets.

Use Instructions

Instructions COMMON to Tandem and Ovoid (Fletcher and Henschke) and Tandem and Ring Applicators

1. Each physician is responsible for the type of sterile technique and type of procedure room utilized for the placement of the applicators.
  2. Perform a thorough pelvic exam.
  3. Place a Foley catheter in the bladder, under sterile conditions per physician discretion.
  4. Determine depth of the uterus by sounding the uterus and taking a measured length.
  5. Select the tandem with the appropriate length and angle based on the patient anatomy and desired intrauterine curvature. The tandem curvature is chosen for each clinical situation.
  6. Optional Smit Sleeves may be used, but not required.
  7. Place the cervical stop at the appropriately measured position on the tandem tip and tighten the hex screw with the hex driver provided. Do not overtighten (hand tighten only).
- Note:** A cervical stop is not used for the Tandem and Ring applicator.
8. It is recommended to assemble the device completely outside of the patient first to assess desired geometry. Make note of the relationship of the attachment ends of the tandem, colpostats, and ring so once the device is inside the patient, the desired geometry can be assured.
  9. Insert the tandem into the uterus to the desired length. The cervical stop should abut the cervix.
  10. Transabdominal ultrasound may be used at physician discretion to verify placement in the uterus.

Instructions UNIQUE to Tandem and Ovoid Applicators (Fletcher and Henschke)

1. Follow common instructions 1-10 above.
  2. The Fletcher applicator is typically used for patients with average anatomy and geometry. The Henschke is typically used for patients with smaller anatomy.
- Note:** For the Fletcher-type applicator, 2 types of ovoids are available, "barrel" shaped (rounded) and "Fletcher" (cylindrical).
3. Select the ovoids of the desired size and attach them to the ends of both the right and left colpostats. Assure the ovoids are attached firmly (hand tighten only) at the distal end of both colpostats.
- Note:** The ovoids for the **Fletcher** set screw onto the end of the colpostats (hand tighten only).
- Note:** The ovoids for the **Henschke** set slide onto the end and are secured with a set screw using the hex driver. The ovoids are notched for proper positioning on the colpostats. Push the ovoid onto the colpostat until it custom fits on the colpostat locating pin. The flat sides of the ovoids face medially.
4. Insert the right colpostat/ovoid (**patient's right side**) into position in the vaginal fornix first, and clip the tandem into the tandem clamp. Secure the locking knob on the underside of the clamp to hold the tandem in place.
  5. Insert the left colpostat/ovoid (patient's left) into position and the proximal end of the spacing bracket will align with the right side and is secured by the locking screws.
  6. When assembling, the line on the tandem clamp should be physician facing.
  7. Tighten both the tandem lock knob (underneath) and the spacing bracket lock knob (on top) to

- ensure the tandem and ovoid assemblies are secured in place in the desired treatment geometry.
8. Vaginal packing (or equivalent) can be placed anteriorly and posteriorly to the device to displace the bladder anteriorly and the rectum posteriorly (as clinically indicated), and secure the device.
  9. Optional: To take radiographs, remove the sterilization caps from the colpostats and tandem and insert the X-ray marker fully into each channel.
  10. Utilize treatment planning images (i.e., CT) according to standard guidelines as required for dosimetry. When the images are completed, remove the X-ray markers (if applicable) and replace the sterilization caps until ready to attach the transfer tubes.

Instructions UNIQUE to Tandem and Ring Applicators

1. Follow common instructions 1-10 above.
  2. Select the ring that corresponds to the degree of anteversion or retroversion of the uterine body and cervical angle.
- Note:** The ring is angled to complement the angle of the tandems. Each angle is clearly indicated on the proximal end of the tube. **Each pair (ring and tandem) is color-coded to match.** The sum of the angles of the rings and tandems equal 90° (i.e., 45+45, 30+60, 60+30).
- Note:** Ensure that the ring angle chosen complements that of the tandem being used (i.e., 45+45, 30+60, 60+30). If incompatible components are used, the applicator will be misaligned and the correct treatment will not be delivered.
3. Fit the ring cap onto the end of the ring and ensure that the build-up cap has the correct orientation (facing the cervix). The two ring cap pieces are screwed together (hand tighten only).
  4. Slide the ring bracket over the proximal ends of the ring shaft until it stops at the colored mark and secure it in place with the upper screw. The closed channel corresponds to the ring tube.
- Note:** The bracket is symmetrical and can be used with the locking screws on the physician's right or left side depending on preference.
5. Slide the ring over the proximal end of the tandem into position against the cervix. The ring cap is beveled, so that it can be angled allowing a smaller profile while entering the patient.
  6. The open, slotted end of the bracket is then lowered onto the tandem and secured in place with the locking screw.
- Note:** If in proper position, the circular markings on the tandem and ring shafts will align for easy verification that the device is assembled correctly.
7. Vaginal packing (or equivalent) can be placed anteriorly and posteriorly to the device to displace the bladder anteriorly and the rectum posteriorly (as clinically indicated), and secure the device.
  8. Optional: To take radiographs, remove the sterilization caps from the tandem and ring and insert the X-ray marker fully into each channel.
  9. Utilize treatment planning images (i.e., CT) according to standard guidelines as required for dosimetry. When the images are completed, remove the X-ray markers (if applicable) and replace the sterilization caps until ready to attach the transfer tubes.

Instructions COMMON to Tandem and Ovoid (Fletcher and Henschke) and Tandem and Ring Applicators

1. Remove the sterilization caps from the tandem and colpostats/ring prior to radiation delivery and attach the transfer tubes from the HDR afterloading device, ensuring the connection between the tandem and colpostats/ring is secure.
2. Ensure there is no tension or excessive bending in the transfer tube.
3. Before starting the treatment, as necessary, inform the patient about appropriate behavior during the treatment. The patient should remain still in the supine position.
4. When treatment is completed, disconnect the transfer tubes from the tandem and colpostats/ring and replace the sterilization caps on the tandem and colpostats/ring.

Removal and Disassembly

1. Gently remove the vaginal packing (if applicable).
2. Remove each applicator in reverse order from placement.
3. Remove the Foley catheter if utilized.
4. Initiate cleaning protocol immediately upon removal of device.

Assembly and Use

*Kobold Vaginal Cylinder Applicator Set*

Assembly Instructions

- The sterilization cap should be in place on the obturator during cleaning and sterilization to prevent foreign materials from entering the tube. It is recommended that the sterilization cap not be removed until the applicator is positioned in the patient and removal is necessary to connect to the HDR afterloader device.
1. Select the vaginal cylinder with the desired diameter based on clinical measurements.
  2. Hand tighten the vaginal cylinder obturator collet into the cylinder chosen for the procedure. (Do not overtighten; hand tighten only.)
  3. Loosely thread the vaginal cylinder obturator collet cap over the collet.
  4. Place the obturator into the central channel and push it up to the tip of the vaginal cylinder until it stops. Verify that the obturator is fully inserted to the tip of the cylinder.
  5. With the obturator held in place, the collet cap can be tightened to lock the obturator in place.
- Warning:** Screwdriver placed through hole of collet should not be used to tighten collet into obturator. Screwdriver may only be used to loosen. Any misuse will void warranty.

Use Instructions

1. Select the vaginal cylinder with the desired diameter based on clinical measurements.
2. Follow Assembly Instructions above.
3. Place the rounded tip of the vaginal cylinder into the vagina with sterile surgical lubrication on the shaft of the cylinder to reduce vaginal friction.
4. The method of holding the cylinder in place is according to the physician's preference; however, securing the device location in the vagina is strongly recommended.
5. Optional: To take radiographs, remove the sterilization cap from the obturator and insert the X-ray marker fully into each channel.
6. Utilize treatment planning images (i.e., CT) according to standard guidelines as required for

dosimetry. When the images are completed, remove the X-ray markers (if applicable) and replace the sterilization caps until ready to attach the transfer tubes.

7. Remove the sterilization cap from the obturator end prior to radiation delivery and attach the transfer tube from the HDR afterloader device, ensuring the connection between the obturator and transfer tube is secure.

8. Ensure there is no tension or excessive bending in the transfer tube.

9. Before starting the treatment, as necessary, inform the patient about appropriate behavior during the treatment. The patient should remain still in the supine position.

10. When treatment is completed, disconnect the transfer tube from the obturator and replace the sterilization cap on the obturator.

#### Removal and Disassembly

1. Gently remove the applicator from the vagina.
2. Loosen the collet cap.
3. Withdraw the obturator from the cylinder.
4. Unscrew the collet cap and collet, removing from the cylinder.
5. Initiate cleaning protocol immediately upon removal of device.

#### Instructions COMMON to All Applicators

*Kobold Fletcher Tandem and Ovoid Applicator Set*

*Kobold Henschke Tandem and Ovoid Applicator Set*

*Kobold Tandem and Ring Applicator Set*

*Kobold Vaginal Cylinder Applicator Set*

#### Storage and Repeated Use

Store the applicator set so it will not be damaged. The applicator set is to be stored in its Kobold-supplied sterilization tray. Ensure the sterilization caps are fitted to the tandem/obturator and ring/colpostats at all times, except when the transfer tube is attached. Store at a temperature between 0° and 30°C.

All parts of the applicator set were designed for repeated use, except the sterilization caps which are single-use only.

**Note:** If the proper condition of a part is in question, that part must not be used. The part should be returned immediately for inspection by the manufacturer. If damage to a part is obvious, please contact the manufacturer for replacement part or servicing.

The manufacturer is not responsible for the use of damaged parts or the use of parts that do not meet factory specifications.

#### Disposal

It is the user's responsibility to dispose of all parts according to national regulations. Upon disposal of the product or related parts, it is suggested to:

- Check the parts for radioactive contamination.
- Disinfect the parts before disposal.
- Dispose sharp objects in a safe way.

#### Warranty

The warranty period for Kobold applicators is for 66 months from the date of completed delivery or 60 months from product acceptance (whichever comes sooner).

Should an applicator fail to conform or otherwise be found to be defective (excluding normal wear and tear) within the applicable warranty period, Kobold's warranty remedy is at Kobold's option repair or replacement of applicator.

Warranty is not valid if equipment is damaged due to customer's mishandling or misuse of said equipment.

Kobold makes no warranty, whatsoever, to any persons other than the original customer.

#### Kobold Brachytherapy Technical Support

For technical questions relating to this applicator, its use, or for reporting any incidents involving the applicator, contact Kobold Brachytherapy Technical Support: 509-703-5090 | [www.koboldmedical.com](http://www.koboldmedical.com) | [info@koboldmedical.com](mailto:info@koboldmedical.com).

#### Source Positions

| WALL THICKNESS (DISTAL END)                           |           |          |
|-------------------------------------------------------|-----------|----------|
| Kobold Applicator Set                                 | Colpostat | Tandem   |
| Kobold Fletcher-model Tandem and Ovoid Applicator Set | 0.396 mm  | 0.396 mm |
| Kobold Henschke-model Tandem and Ovoid Applicator Set | 1.9812 mm | 0.396 mm |
|                                                       | Ring      | Tandem   |
| Kobold Tandem and Ring Applicator Set                 | 0.396 mm  | 0.396 mm |
| Obturator                                             |           |          |
| Kobold Vaginal Cylinder Applicator Set                | 0.6662 mm |          |

Manufactured For:  
Kobold, LLC  
23403 E. Mission Ave., Suite 220E  
Liberty Lake, WA 99019  
USA

Authorized Representative:  
Emergo Europe  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands

United States Patent No. 9,180,311  
United States Patent No. 9,895,551

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